

Supplementary Table 1. Primers used for genotyping and PCR analysis

	Forward (5'-3')	Reverse (5'-3')
Genotyping		
SOCS2	TGTTTGACTGAGCTCGCGC	CAACTTTAGTGTCTTGGATCT
Neo	ACCCTGCACACTCTCGTTTTG	CCTCGACTAAACACATGTAAAGC
PCR Analysis		
<i>Socs1</i>	TCCGATTACCGGCGCATCACG	CTCCAGCAGCTCGAAAAGGCA
<i>Socs2</i>	TGGCTGCTCAAGATCAAATG	TGTCCTCCTGGAAATGGAAG
<i>Socs3</i>	GAGTACCCCAAGAGAGCTTACTA	CTCCTTAAAGTGGAGCATCATACTG